



RENPRO 188R

Pro-Grade Athletic Tape

DESCRIPTION

The Pro-Grade Athletic Tape provides support and helps to restrict motion in muscles and joints. This tape is made from 100% cotton and has a zinc-oxide (skin-friendly), vanilla-scented adhesive. Pair it up with our Foam Under Wrap for extra cushioning.

APPLICATIONS

- Wrists/Ankles

PRODUCT BENEFITS

- Helps to prevent injury or support existing injuries
- Pairs well with Renfrew Foam Under Wrap
- Resists failure caused by moisture

TECHNICAL PROPERTIES

Technical Property	Nominal Value	Unit	Test Method
Total Thickness	8.5 (.215)	Mils (mm)	REN TM 20
Adhesion to Steel	28 (7.7)	Oz/Inch (N/25mm)	PSTC - 101 Mod
Adhesion to Backing	16 (4.4)	Oz/Inch (N/25mm)	PSTC - 101 Mod
Tensile Strength	34 (151)	lbs. / inch (N/25mm)	ASTM D 1000

This Technical Data Sheet (TDS) is provided by Scapa for information only. The contents of this TDS are considered to be reliable but accuracy is not guaranteed. The results shown in this TDS may not be relevant to all users of the product and it is the user's responsibility to determine the suitability of the product for its intended purpose. Scapa expressly disclaims all express or implied warranties including merchantability and fitness for a particular purpose to the fullest extent permitted by law and does not accept responsibility or liability (whether direct, indirect, or consequential) arising from use of the product or reliance on the TDS. Scapa does not exclude liability for death or personal injury arising from its negligence or fraud. Scapa refers to Scapa Group plc and its subsidiaries.

STANDARD PRESENTATIONS

- Colours: White

RECOMMENDATIONS

Application:

Renfrew Athletics trainers tape may be used in conjunction with other Renfrew Athletics products to provide a custom applied stabilization while maintaining complete comfort to the athlete. In particular, Renfrew Athletics recommends pre-wrapping any area to be taped with our foam underwrap tape for maximum comfort.

Athletic tapes should only be applied by trained professional. Incorrect wrapping may not provide the support needed by the athlete, and if applied too tight, may restrict circulation.

NOTE: Shelf life is one year from date of shipment, when stored in a cool dry place below 76°F (24°C). Rolls should be stored on end.